


EXPORTING FROM **CARESTREAM**

UPLOAD GUIDE



866-497-3692 • glidewelldental.com

Note: If using the Instrumentarium Dental OP300 or Gendex GXDP-700™, instructions will not apply. Please contact your local representative for assistance on extracting the raw data from your CT scanner.



*Step 1: Right-click on an empty space on your desktop to create a new folder (New → Folder).
Type in the patient name and click “Enter”.*

1




Step 2: Open image and go to “Export the volume” icon.

Note: If using the Instrumentarium Dental OP300 or Gendex GXDP-700™, instructions will not apply. Please contact your local representative for assistance on extracting the raw data from your CT scanner.



Step 3: In the pop-up window for **1. Select file and directory**, click “Browse” to select the export location.

2



Step 4: Select your newly created folder on the desktop and click “OK”.

Note: If using the Instrumentarium Dental OP300 or Gendex GXDP-700™, instructions will not apply. Please contact your local representative for assistance on extracting the raw data from your CT scanner.

1. Select file and directory

Select output directory: C:\Users\D15189\Desktop\Test, Test [Browse]

2. Select resolution

Select axial matrix size (pixels):
☒ 543x543 - Same as source
☐ 512x512
☐ 256x256

Select slice thickness:
☐ Create isotropic volume
 0.150 mm - Same as source


3. Select output file format

Select file format:
☒ CT (one file per slice)
☐ Multi-frame (one single file)
☐ Anonymize
☐ Create DICOMDir

Select compression:
☒ Uncompressed
☐ Compressed (JPEG lossless)


Export

Step 5: No change for **2. Select resolution**, and for **3. Select output file format**, click on “CT (one file per slice)” and “Uncompressed” and then “Export”.



Step 6: Go to the patient's folder and confirm the .dcm files were exported correctly (there should be several hundred files).

Note: If using the Instrumentarium Dental OP300 or Gendex GXDP-700™, instructions will not apply. Please contact your local representative for assistance on extracting the raw data from your CT scanner.



Step 7: Return to your desktop and right-click on the patient's folder to send to a compressed folder (Send to → Compressed (zipped) folder). A duplicate folder will be created with either a zipper or a blue "Z". Note: If following a Dual Scan Protocol, zip scans individually and label accordingly (i.e., "Patient Scan" and "Dental Scan").

4




Step 8: Upload your case on glidewelldental.com, then click and log in to "My Account". Note: Use Google Chrome when uploading your case; additionally, if you do not have an account, please contact us and we will set it up for you. Your scan can also be sent on a disc or USB drive via standard shipping methods if digital upload is not possible.

Note: If using the Instrumentarium Dental OP300 or Gendex GXDP-700™, instructions will not apply. Please contact your local representative for assistance on extracting the raw data from your CT scanner.




Step 9: Select "Cases" → "Guided Surgery/DTP".



Step 10: Select "Request New Plan".

5



Step 11: Select desired package.