


NEW!

Beauty

in a Class of Its Own



Here's your guide to BruxZir® Esthetic, the newest solution from the **No.1** prescribed brand of full-contour zirconia

BruxZir®
Esthetic
—SOLID ZIRCONIA—

NEW!

BruxZir[®]
Esthetic
— SOLID ZIRCONIA —

Beauty
in a Class of Its Own

- ▶ Optimized combination of strength and translucency
- ▶ Lifelike vitality, engineered for the anterior
- ▶ Enhanced shade matching for predictable results

BruxZir[®] Esthetic Solid Zirconia is a new category of dental zirconia made possible by years of research & development at Glidewell Dental. **Named in Dentistry Today's Top 100 Products in 2018**, BruxZir Esthetic zirconia is a revelation for anterior cases and further proof of why BruxZir Solid Zirconia is the **No. 1 prescribed brand of full-contour zirconia**.



See back cover for more product details

DELIVER STUNNING
RESULTS WITH
BRUXZIR[®] ESTHETIC




This patient's severe anterior wear and loss of canine and protrusive guidance were chief concerns. To provide strong, lifelike restorative results, BruxZir Esthetic zirconia was chosen for a 3-unit bridge for teeth #11-13, crowns for teeth #21 & #22, and veneers for teeth #14 & #23.

— Clinical dentistry by Dr. Anamaria Muresan

TRUE-TO-NATURE TRANSLUCENCY


BruxZir Esthetic Outperforms Other Zirconias



The translucency of BruxZir Esthetic zirconia is superior to that of other zirconias. High translucency in the warm color wavelengths enables **more lifelike restorations** with BruxZir Esthetic zirconia.

HIGH-PERFORMING STRENGTH


BruxZir Esthetic Sets a New Level of Durability




At 870 MPa, BruxZir Esthetic zirconia provides an average flexural strength that is **74 percent stronger than lithium disilicate** and **over 33 percent stronger than BruxZir Anterior zirconia**.

SUPERIOR PROCESSING FOR **ULTRA PERFORMANCE**

SEM of sintered
isostatically pressed zirconia



SEM of sintered colloiddally
processed BruxZir Esthetic



Proprietary processing yields a **smaller, more uniform zirconia grain size**, effectively elevating BruxZir Esthetic zirconia to a class of its own. This specialized processing is key for achieving the material's exquisite translucency and extraordinary strength.

NEW!

BruxZir®
Esthetic
—SOLID ZIRCONIA—

Beauty
in a Class of Its Own

From dentistry's *leading authority* in zirconia technology comes BruxZir® Esthetic Solid Zirconia, a *unique formulation* scientifically engineered for the optimal combination of *strength and esthetics*.

Brilliant Smiles Start with BruxZir® Esthetic


- ▶ Premium-quality monolithic zirconia
- ▶ High-performing flexural strength of 870 MPa (average)
- ▶ Lifelike translucency
- ▶ New, superior shading technology for enhanced shade matching and consistency
- ▶ Ideal for anterior crowns & bridges of up to 4 units
- ▶ Advanced, proprietary colloidal processing



All Restorations
Made in México

Preparation Guidelines

BruxZir Esthetic Solid Zirconia



For more information
mx.glidewell dental.com ■ 01-800-212-9080

 **GLIDEWELL
LABORATORIOS**

MKT 11988_1.0 GLMX-757200-031119 EN